

30° Sliding Binocular Head

Type: Sliding Binocular Head

Viewing Angle: 30° inclined 360° Rotatable

Eyeiece Tubes: Two, for binocular vision

Interpupillary Distance: Adjustable (typically 48-75 mm)

Diopter Adjustment: Individual for each eyepiece

Sliding Mechanism: Allows horizontal adjustment of eyepiece tubes for customized positioning

Magnification: Compatible with standard eyepieces (e.g., 10x, 15x, 20x)

Optics: High-quality optics for clear, bright images

Material: Durable metal or high-grade plastic

Mounting: Fits microscope mounts suitable for 40mm or 42mm Dia.
(Any One as per requirement)

This head offers flexible viewing and ergonomic adjustments for comfort during extended use.

